





SonoTip® EBUS Olympus and Fujifilm

The Right Needle System for Every Situation



When Results Count:

- NEW: Compatibility with all Olympus bronchoscope types
- For obtaining high quality biopsies
- High visibility under ultrasound



SonoTip® EBUS Pro & EBUS ProFlex

TBNA with Facet Cut

For atraumatic punctures low resistance and large needle opening for adequate material uptake

Compatible with Olympus (type BF-UC180F and BF-UC190F) / Fujifilm

With atraumatic (round) nitinol stylet; click-lock adapter, incl. VacLok™

Aspiration Syringe with a filling volume of 5 - 20 ml and Olympus BF-UC190F standard adapter; incl. separately packed Olympus BF-UC180F and Fujifilm adapters



Catalogue No.	Description	Needle Material	Needle Diameter (Gauge/mm)	Puncture Depth of Needle (cm)	Outer Tube Diameter (cm)
GUS-45-18-022	SonoTip® EBUS Pro: EBUS guided TBNA system	Stainless Steel	22/0,7	0-4	1,8
GUS-46-18-022	SonoTip® EBUS ProFlex: EBUS guided TBNA system	Nitinol	22/0,7	0-4	1,8

SonoTip® TopGain EBUS

TBNB with 3-Point Crown-Cut

Dedicated new needle design for tissue biopsy -
Allows the extraction of coherent tissue to support the histological evaluation

Compatible with Olympus (type BF-UC180F and BF-UC190F) / Fujifilm

With atraumatic (round) nitinol stylet; click-lock adapter, incl. VacLok™


Aspiration Syringe with a filling volume of 5 - 20 ml and Olympus BF-UC190F standard adapter; incl. separately packed Olympus BF-UC180F and Fujifilm adapters



Catalogue No.	Description	Needle Material	Needle Diameter (Gauge/mm)	Puncture Depth of Needle (cm)	Outer Tube Diameter (cm)
GUB-45-18-022	SonoTip® TopGain EBUS: EBUS guided TBNB system	Stainless Steel	22/0,7	0-4	1,8
GUB-46-18-022	SonoTip® TopGain EBUS: EBUS guided TBNB system	Nitinol	22/0,7	0-4	1,8



SonoTip® Needle Features




Needle in Two Variants – Stainless Steel and Nitinol

Reliable stainless steel needle with proven puncture behavior and 100% Nitinol needle with Shape Memory Effect for improved flexibility

NEW: Combination Adapter - Suitable for All Olympus Bronchoscopes

- Compatibility with Olympus bronchoscope types **BF-UC190F & BF-UC180F**
- Ensures stable connection to the working channel of the bronchoscope - no wobbling, no slipping
- With „click“ sound to confirm secure fixation





Visit our website:
www.medi-globe.de/en



Medi-Globe GmbH
on LinkedIn



SonoTip® EBUS Olympus & Fujifilm
on Youtube

Olympus BF-UC190F and
BF-UC180F Adapter



Fujifilm-Adapter



Medi-Globe GmbH

Medi-Globe-Straße 1–5 • D-83101 Rohrdorf OT Achenmühle
Tel.: +49 (0) 8032 973-379 • Fax: +49 (0) 8032 973-392
E-Mail: sales@medi-globe.de • www.medi-globe.de